

Non-invasive, drug free
therapy for Sports



SCENAR *for Sports*

www.scenar.com.au

Self Controlled Energo Neuro Adaptive Regulator

Designed for Professional athletes

Sports people

Coaches

Fitness instructors



Try SCENAR SPORT - Beat the competition

RITMSCENAR Sport is an electro-therapy device with built-in electrode that delivers complex electrical signals through the skin. These signals stimulate the brain to generate Neuropeptides, the key biochemicals needed by the body to heal itself.


Neuropeptides may be considered as our body's own internal pharmacy.

RITMSCENAR Sport uses a sophisticated software program which controls and automatically changes the signal waveforms and frequencies to prevent the body from adapting to the stimulation.

RITMSCENAR Sport is designed for sport professionals to increase the sportsmen effectiveness, reduce the trauma downtime and fasten the rehabilitation process.

Pain, cramping soft tissue injuries, contusions, fractures, tendon pulls and ligament tears are common conditions readily treatable by SCENAR therapy. Treatment after event has also proven to be effective in recovering from tired and aching muscles. Physical exhaustion, and other post-event aches and pains.

Using **RITMSCENAR Sport** does not require medical background, so every coach and sportsman can benefit from this exceptional device derived from the Russian Space technology.



“SCENAR supported my nervous system, encouraged my body to heal, allowed me to maintain strength and fitness so I could compete at my optimum level ”



Layne Beachley
7 times World Surfing Champion
Australia

“The SCENAR technology is the new sensation in sports performance.

We have used it during our preparation for the 2008 Rowing Season and the Junior World Championships.

The SCENAR technology was very effective in speeding the recovery process, injury prevention and quick pain relief. The feedback from the athletes was positive and their general health improved.

SCENAR is incredible technology with wide application in Sports. ”



Stani Slavova
NSW Elite Development Coordinator – Rowing
Coach of the 2008 Australian Rowing Junior Team

Certificates:



Australia
ARTG
ID: 140659

U.S.A
FDA
882.5050

Device Parameters

Diagnostic Modes	Integral dose
Dose indication	Sound and light
Amplitude modulation	3 : 1
Frequency modulation	30 -120 Hz, cycle 7 sec
Influencing impulse frequencies	14, 60, 90, 340 Hz
Intensity <i>(the number of impulses in a packet)</i>	1
Amplitude of the stimulating pulse at a standard load	(1,7-2,8) ... (100-140) V
Indicator	Light emitting diodes
Power supply	one 9V alkaline battery
Temperature	10 to 35C (Operating) -50 to 50C (Storage & Transportation)
Weight	0.2kg
Dimensions	130 x 50 x 25mm
Warranty	2 years



Manufacturer – Head Office
RITM OKB ZAO

Petrovskaya 99,
Taganrog 347900
Russian Federation

Tel/fax: +7 8634 623179
www.scenar.com.ru
medsc@scenar.com.ru

TRAINING

Specialized training for SCENAR Sport users is available. It covers theoretical and practical applications of the device and technology.

For more information about training options, please contact the Authorized Training Centre:



Australia, New Zealand, Oceania
RITM Australia Pty Ltd

Unit 4,
130-134 Pacific Hwy
Greenwich NSW 2065
Australia

Tel: +61 2 8011 4217
Fax: +61 2 9439 1930
info@scenar.com.au
www.scenar.com.au



Global Scenar Academy Pty Ltd.
7/12 Rickard Road
Narrabeen NSW 2101 Australia
Tel: +61 2 9913 2020
info@scenartraining.com
www.scenartraining.com

Distributor: _____

Self Controlled Energo Neuro Adaptive Regulator

SCENAR *for Sports*